

ST-LED60S

ST-LED60M

ST-LED60E

TECHNICAL DATA

Model No.	ST-LED 60D	ST-LED 60S	ST-LED 60M	ST-LED 60E
Bulb	LED	LED	LED	LED
Diameter of Lamp Head	60 cm + 60 cm	60 cm	60 cm	60 cm
Light Intensity in Lux at 1 Meter	140,000 Lux + 140,000 Lux	140,000 Lux	140,000 Lux	140,000 Lux
Color Temp.	4200 ± 200°K	4200 ± 200°K	4200 ± 200°K	4200 ± 200°K
Diameter of Illuminated Field at 1 M	20 cm + 20 cm	20 cm	20 cm	20 cm
Color Rendition Index	90	90	90	90
Working Distance	60 – 150 cm	60 – 150 cm	60 – 150 cm	60 – 150 cm
Bulb Life / Hours	> 30,000	> 30,000	> 30,000	> 30,000
Power Consumption	90W + 90W	90W	90W	90W
Lamp Head High Adjustment	90 cm	90 cm	130 cm	130 cm
Electrical Connection Input / AC	110 – 230V	110 – 230V	110 – 230V	110 – 230V
Minimum Ceiling High (An Extension Column is Available For Higher Ceiling at Extra Cost)	240 cm	230 cm	---	---
Total Weight / Approx	120 Kgs	60 Kgs	110 Kgs	80 Kgs

ST-LED 60

NEW LIGHT, NEW CHOICE



ST. FRANCIS Medical Equipment Co., Ltd.
TAIPEI TAIWAN R. O. C.

CE ISO9001 ISO13485 FDA, 510(k)
GMP CERTIFIED FACTORY