

**COMPARATIVE CHART OF OT-2300, OT-2200, OT-2100 ,OT-2000 AND OT-2000L (2019JAN.Update)**

|                                                                | <b>OT-2300</b>                                                                    | <b>OT-2200</b>                                                                      | <b>OT-2100</b>                                                                      | <b>OT-2000</b>                                                                      | <b>OT-2000L</b>                                                                     |
|----------------------------------------------------------------|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
|                                                                |  |  |  |  |  |
| <b>SPECIFICATIONS</b>                                          |                                                                                   |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
| <b>Table Dimension</b>                                         |                                                                                   |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
| Length                                                         | 2125 mm                                                                           | 2160 mm                                                                             | 2160 mm                                                                             | 2160 mm                                                                             | 2160 mm                                                                             |
| Width                                                          | 600 mm                                                                            | 600 mm                                                                              | 600 mm                                                                              | 600 mm                                                                              | 600 mm                                                                              |
| Table Top Elevation(with mettress) / Electric Hydraulic        | 725~1070 mm                                                                       | 700~1000 mm                                                                         | 740~1040 mm                                                                         | 740~1040 mm                                                                         | 560~970 mm                                                                          |
| Trendelenburg/ Reverse Trendelenburg / Electric Hydraulic      | 30°/ 30°                                                                          | 20°/ 20°                                                                            | 20°/ 20°                                                                            | 30°/ 30°                                                                            | 25°/ 25°                                                                            |
| Lateral tilt (Left/ Right) / Electric Hydraulic                | 30°/30°                                                                           | 20°/20°                                                                             | 20°/20°                                                                             | 20°/20°                                                                             | 20°/20°                                                                             |
| Lower Back Section (Up/Down) / Electric Hydraulic              | 85°/50°                                                                           | 80°/30°                                                                             | 80°/30°                                                                             | 80°/45°                                                                             | 80°/45°                                                                             |
| Upper Back Section (Up/Down) / Electric Hydraulic              | 65°/60°                                                                           | N.A                                                                                 | N.A                                                                                 | N.A                                                                                 | N.A                                                                                 |
| Head Section (Up/Down) Manual                                  | 25°/90°                                                                           | 25°/90°                                                                             | 25°/90°                                                                             | 25°/90°                                                                             | 25°/90°                                                                             |
| Foot Plate (Up/Down)                                           | 0°/90° (Electric Hydraulic)                                                       | 20°/90°(Manual)                                                                     | 20°/90°(Manual)                                                                     | 20°/90°(Manual)                                                                     | 20°/90°(Manual)                                                                     |
| Flex/ Reverse Flex / Electric Hydraulic                        | 205°/80°( <b>Auto-setting</b> )                                                   | 210°/90° (Auto-setting)                                                             | N.A                                                                                 | N.A                                                                                 | N.A                                                                                 |
| Split Angle / (Foot Plate) / Manual                            | N.A                                                                               | 245°                                                                                | 245°                                                                                | 245°                                                                                | 245°                                                                                |
| Table Top Sidewise Slide (Left/Right) / Electric Hydraulic     | N.A                                                                               | 150 mm                                                                              | N.A                                                                                 | N.A                                                                                 | N.A                                                                                 |
| Table top longitudinal slide (Left/Right) / Electric Hydraulic | N.A                                                                               | 350 mm                                                                              | 350 mm                                                                              | N.A                                                                                 | N.A                                                                                 |
| 4 electro-hydraulic floor locking castors                      | YES                                                                               | YES                                                                                 | YES                                                                                 | YES                                                                                 | YES                                                                                 |
| Auto-Resetting "0" position of table top                       | YES                                                                               | YES                                                                                 | YES                                                                                 | YES                                                                                 | YES                                                                                 |
| <b>ACCESSORIES</b>                                             |                                                                                   |                                                                                     |                                                                                     |                                                                                     |                                                                                     |
| Emergency Override Panel                                       | Standard                                                                          | Standard                                                                            | Optional                                                                            | Optional                                                                            | Optional                                                                            |
| UPS Battery Device                                             | Standard                                                                          | Standard                                                                            | Optional                                                                            | Optional                                                                            | Optional                                                                            |
| Hand controller                                                | Standard                                                                          | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Armrest/2 pcs                                                  | Standard                                                                          | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Screen Frame                                                   | Standard                                                                          | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Knee crutches                                                  | Option                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Headrest                                                       | Standard                                                                          | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Shoulder support/2 pcs                                         | Option                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Body Support                                                   | Option                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Mattress                                                       | Standard                                                                          | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |
| Lateral support                                                | Option                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            | Standard                                                                            |

|                                                                            |              |              |              |              |              |
|----------------------------------------------------------------------------|--------------|--------------|--------------|--------------|--------------|
| <b>Additional Accessories ( option )</b>                                   |              |              |              |              |              |
| Wireless remote Controller(no. Ot-1832)                                    | YES          | N.A          | N.A          | N.A          | N.A          |
| Foot support (no. Ot-1724)                                                 | YES          | YES          | YES          | YES          | YES          |
| Cylindrical support/s.s. (no. Ot-1726)                                     | YES          | YES          | YES          | YES          | YES          |
| Raised armrest/s.s. (no. Ot-1728)                                          | YES          | YES          | YES          | YES          | YES          |
| Funnel and wast receptacle and corrugated tubing/s.s. (no. Ot 1718)        | YES          | YES          | YES          | YES          | YES          |
| U-shaped headrest/s.s. (no. Ot-1720)                                       | YES          | YES          | YES          | YES          | YES          |
| Head cushion (no. Ot-1723)                                                 | YES          | YES          | YES          | YES          | YES          |
| Kidney elevator (no. Ot-1740)                                              | YES          | YES          | YES          | YES(OT-1624) | YES(OT-1624) |
| Waterproof Foot Pedal controller                                           | YES          | N.A          | N.A          | N.A          | N.A          |
| Wrist support (For Ophthalmic position)                                    | YES          | YES          | YES          | YES          | YES          |
| Sitting head support (no Ot-1732) (For Sitting position)                   | YES          | YES          | YES          | YES          | YES          |
| Chest/abdominal support (no. Ot-1730)                                      | YES          | YES          | YES          | YES          | YES          |
| Instrument tray/s.s. (no. Ot-1725)                                         | YES          | YES          | YES          | YES          | YES          |
| Abdominal support (no. Ot-1727)                                            | YES          | YES          | YES          | YES          | YES          |
| Table top x-ray radiography sheets (no. Ot-1729)                           | YES          | YES          | YES          | YES          | YES          |
| Table for hand and arm surgery/s.s. (no. Ot-1719)                          | YES          | YES          | YES          | YES          | YES          |
| Restraint strap (no. Ot-1620)                                              | YES          | YES          | YES          | YES          | YES          |
| Infusion rod/s.s. (no. Ot-1721)                                            | YES          | YES          | YES          | YES          | YES          |
| Orthopaedic device with cart for lower limbs/s.s. (no. Ot-090)             | N.A          | YES          | YES          | YES          | YES          |
| Orthopaedic device with cart for lower limbs and forearm/s.s. (no. Ot-080) | N.A          | YES          | YES          | YES          | YES          |
|                                                                            |              |              |              |              |              |
| <b>Physical Properties</b>                                                 |              |              |              |              |              |
| Table weight (approx) kgs                                                  | 285          | 285          | 285          | 170          | 285          |
| Mattress Thick                                                             | 5.5 cm       | 5.5 cm       | 5.5cm        | 5.5 cm       | 5.5 cm       |
| Patient Weight capacity                                                    | 250 Kgs.     | 250 Kgs.     | 250 kgs.     | 180 kgs.     | 250 kgs.     |
|                                                                            |              |              |              |              |              |
| <b>Power requirements</b>                                                  |              |              |              |              |              |
| Electrical AC Input                                                        | 110~230V     | 110~230V     | 110~230V     | 110~230V     | 110~230V     |
| Frequency                                                                  | 50/60 Hz     | 50/60 Hz     | 50/60 Hz     | 50/60 Hz     | 50/60 Hz     |
| Remote Control Power(no. Ot-1832)                                          | 12 Vdc       | N.A          | N.A          | N.A          | N.A          |
| Safety Standards                                                           | CE, FDA, GMP | CE, FDA, GMP | CE, FDA, GMP | CE, FDA, GMP | CE, FDA, GMP |